



TABELLA MISURE SILOSHEAT

SILOSHEATH TABELLA MISURE						
Mis.	Code	Description	Approx. Gel Length	Approx. Gel Length	Approx. Distal Circ.	Approx. Prox. Circ.
						
I°	13315	Child	20 cm	25 cm - 34 cm	17 cm - 24 cm	17 cm - 26 cm
II°	13325	Short / Narrow	25 cm	38 cm - 44 cm	18 cm - 28 cm	18 cm - 32 cm
III°	13335	Short / Regular	25 cm	38 cm - 44 cm	20 cm - 34 cm	20 cm - 38 cm
IV°	13345	Medium / Narrow	33 cm	43 cm - 49 cm	20 cm - 34 cm	20 cm - 38 cm
V°	13355	Medium / Regular	25 cm	43 cm - 49 cm	23 cm - 37 cm	23 cm - 43 cm
VI°	13365	Large / Narrow	33 cm	46 cm - 54 cm	23 cm - 37 cm	23 cm - 43 cm
VII°	13375	Large / Regular	33 cm	46 cm - 54 cm	28 cm - 42 cm	28 cm - 49 cm
VIII°	13385	Large / Wide	33 cm	46 cm - 54 cm	31 cm - 50 cm	31 cm - 56 cm



COME INDOSSARE IL COPRIMONCONE PROTESICO CON GEL SILOSHEATH



COME TOGLIERE IL COPRIMONCONE